

Nanopore experiment samples

Sample nr.	Description	Type	Methylation enzymes
1	DNA 16*601 (50) to test Hia5	Bare DNA	Hia5
2	DNA 16*601 (50) to test Hia5	Bare DNA	x
3	DNA 16*601 (50) to test Hia5 (with other enzymes)	Bare DNA	M.Cvipi +M.SssI
4	DNA to test Hia5 (with other enzymes)	Bare DNA	Hia5 + M.Cvipi +M.SssI
5	Reconstituted DNA for Hia5 testing (16*601, 50bp)	Reconstitution	Hia5
6	Reconstituted DNA for Hia5 testing (16*601, 20bp)	Reconstitution	Hia5
7	Reconstituted DNA for Hia5 testing (16*601, 50bp)	Reconstitution	Hia5 + M.Cvipi +M.SssI
8	Reconstituted DNA for Hia5 testing (16*601, 20bp)	Reconstitution	Hia5 + M.Cvipi +M.SssI
9	Reconstituted DNA for Hia5 testing (16*601, 50bp) (suboptimal, not enough nucleosomes)	Reconstitution	Hia5 + M.Cvipi +M.SssI
10	Wt yeast grown in GAL	Native chromatin	M.Cvipi +M.SssI
11	Wt yeast grown in RAF	Native chromatin	M.Cvipi +M.SssI
12	Wt yeast grown in RAF	Native chromatin	Hia5

Position —	Flow cell ID —	Sample ID —	Health	Run time —	Run state —	Reads —	Bases —	Basecalled % —
MN30914	ACR467 FA-01378	no_sample	—	1 d 14 h	Complete	8.08 k	34.36 Mb basecalled 31.13 Mb estimated	100%

Scroll right >

Position MN30914

Experiment group 2022-01-18_12samplemethylationtest
Sample ID no_sample
Flow cell product code FLO-FLG001
Kit ID SQK-LSK109
Current output directory D:\Nanopore\2022-01-18_12samplemethylationtes...
Basecall model Modified basecalling for 5mC
Minimum qscore 8
Barcoding EXP-NBD104


 **Total run time** 1d 14h 25m 30s

Read count 8.08 k

Estimated bases 31.13 Mb

Basecalled bases 34.36 Mb

Cumulative output



● Total reads ● Passed reads ● Failed reads ● Skipped reads

Reads Bases



▶ Resume

⏸ Pause

■ Stop

⚙ Start MUX scan

↓ Export PDF


≡

⌵ Experiment group

Position —	Flow cell ID —	Sample ID —	Health	Run time —	Run state —	Reads —	Bases —	Basecalled % —
MN30914	ACR467 FA-01378	no_sample	—	1 d 14 h	Complete	8.08 k	34.36 Mb basecalled 31.13 Mb estimated	100%

Scroll right >

Barcode hits



Sort

Reads

Bases

☐ Display failed☐ Display unclassified☒ Hide zero values

↻ Reset selection